



EBI FOOD SAFETY

### **PCR of phage DNA:**

Phage DNA does not need to be extracted from bacteriophage, rather phages can be used directly in the PCR reaction.

Protocol:

1. Remove a plaque from the top-agar with a Pasteur pipette and transfer to an eppendorf cap with 200  $\mu$ l of SM-phage buffer.
2. Incubate at 30°C for 30 min while shaking.
3. Use 1 $\mu$ l of the supernatant as template for the PCR reaction.

PCR settings:

Cycle 30x  
30 sec 94°C denaturing  
60 sec 50-53°C annealing  
60 sec 68-72°C elongation

Using the following primers:

Forward: 5'-ccttcacgcatcttggttacag (binds P100 genome bp:108867-108888)

Reverse: 5'-caggggtgtatttaggtactc (binds P100 genome bp: 109957-109937)

Samples containing P100 will result in fragments of 1090bp length being generated.

By Dr. Steven Hagens

August 2007

EBI FOOD SAFETY B.V.

NIEUWE KANAAL 7P . 6709 PA WAGENINGEN . THE NETHERLANDS . TEL +31 (0)317 421 414 . FAX +31 (0)317 410 055  
E-MAIL : MAIL@EBIFOODSAFETY.COM . INTERNET: WWW.EBIFOODSAFETY.COM . BANK: 1039-55.844 . IBAN: NL53 RABO 0103 9558 44

EBI FOOD SAFETY IS REGISTERED WITH THE CHAMBER OF COMMERCE HAAGLANDEN UNDER NUMBER 272.79042